Leamington Mennonite Home Long Term Care

INFECTION CONTROL POLICY AND PROCEDURE

CATEGORY: Screening & Immunization	SUBJECT: Vaccine Refrigerator Temperature Monitoring		SECTION: D
DATE:	Administrator's Signature:	J. hr.	POLICY: 13
July 12, 2022 REVISION DATES: November 2023		Lloppola, RN.	

VACCINE REFRIGERATOR TEMPERATURE MONITORING

POLICY:

Refrigerators used for the storage of vaccines will be located in a secure area, accessible to staff only, and will be used only for the storage of vaccines. No personal items are to be placed in any of the refrigerators.

PROCEDURE:

The DNPC will:

- 1) Take refrigerator temperature twice daily. **Refrigerator temperature acceptable range:** 2°C to 8°C
- 2) Record temperature in a logbook supplied by the local Public Health Unit (PHU).
- 3) Only store vaccines supplied by the PHU in the vaccine refrigerator.
- 4) Store vaccines on the middle shelves of the refrigerator, and never on the door or in the drawers.
- 5) Retain the logbook for a period of one year.
- 6) Notify the local PHU if it is noted that the temperature falls outside of the above guidelines. Public Health will advise what is to be done with the vaccine.
- 7) Take immediate steps to repair/replace the refrigerator and re-order replacement vaccine.
- 8) Verify Vaccine expiration dates monthly.
- 9) Dispose of any vaccine noted to be out of date through the local PHU.
- 10) Dispose of any open vial after one month as directed.
- 11) Monitor that only vaccines supplied by the PHU are stored in the vaccine refrigerator.

- 12) Use an insulated container as specified by Public Health to transport vaccines from Public Health to site and while conducting on-site clinics.
- 13) Record temperatures of the insulated container as outlined in referenced guidelines.